



Pulsed Magnetic Field for Neutron and X-ray

H. Nojiri, K. Ohyama and Y. Matsuda
IMR, Tohoku University, Sendai, Japan

High Magnetic Field Spin Science in 100 T -High Field Neutron Experiments




Various spectrometers
Single Crystal
Powder diffraction
Polarization analysis
Inelastic



Why magnetic field?

Couple to Spin and Orbital motion of electron



Precise Control
Soft
Time Structure

Low Field

Control of Spin direction



At 100 T


e 135 K
Mn³⁺ 540 K
Dy³⁺ 1000 K

Phase
control


Non Mag.-Mag.

Structure


MI-transition




Drastic
change




Various High field devices




Miniature Magnet and Portable Capacitor Bank



Small Energy




Small Space





High Magnetic Field

ND in repeating pulsed magnetic field



- Pulse Neutron
- KENS in Tsukuba
- 25 T experiment




Mini Coil for Neutron

The Highest Magnetic Field for ND

Repeating : **25T** @ KEK
(Motokawa et al. 1989)
ND over 20T is still Difficult.

Reactor Facility
High Beam Flux
Beam Focusing is easier.



Development of **Easier and more Diffusive Techniques**


Target : **H=40T**

Goal: Observation of Magnetic Transitions by ~**100 Shots**
Experiments

Development to J-PARC




Double Magnet System




$L_1 \sim 0.1\text{mH}$, $L_2 \sim 1\text{mH}$
 $C = 0.96\text{mF}$ (Short Pulse)
 $C = 3.36\text{mF}$ (Long Pulse)

Longer Pulse
Low loss
Effective injection
Current limit
Design Freedom



Experimental

The Triple Axis Spectrometer, **AKANE** of Tohoku Univ.
@ JRR-3M, JAEA, Japan



$$\lambda = 2.0 \text{ \AA}$$


Colimation: guide-Open-S-Blank-Blank

Scattering Plane: a^*-c^*


H//c




The coil & Sample were cooled in a liquid-He cryostat.






Sample Space
($\phi 5\text{mm}$ L10mm)



Magnet Coil (Cu-Ag)




MnF_2



Conectors

Choke Coil




Compact Experiments!



AKANE

AKANE Data
Acquisition System

Condenser Bank
3.36mF




oscilloscope

Controller

PC

Observation of Spin-Flop at 10T

Integration of (100) magnetic reflection after 100shots




We succeeded in observing intensity change due to SF transition at 10T after 100 shots.

Pulse Interval: ~17min.

100 Pulse ~ 25hr.

Problems to be solved



- Long Duration of experiments: ~25hr
⇒ Higher Beam Flux
- Coil Optimization, low resistive
- Cooling efficiency, Nitrogen coil
- Peak Profiles cannot be observed at the present.
⇒ PSD system




Practical Experiments

Enhancement of Efficiency



Beam Focusing by Neutron Super mirror




~3 times
enhancement in
real space



Position Sensitive Neutron Detector(PSD) System



Profile of 100 of MnF_2




PSD can observe **Peak Profiles**
without step scan.

Short Summary

- We **succeeded** in test experiments in MnF_2 by ND after **100 shots** with a **20T** long pulsed magnet in a reactor facility.
- We also succeeded in generating **30T** pulsed magnetic fields stably and safely on the spectrometer.
- By reducing the resistivity of the coils and Increasing the capacitance of the condenser, We can generate 40T magnetic Fields.

Cryogen-free magnet



18T Cyrogen Free Magnet

18.1T was successfully achieved!



Maximum stress 53MPa@Bi2223 coil
 258MPa@Nb₃Sn middle coil

Specifications

Magnet

- Bi2223 inner most coil (2.7T)
 - co-winding with SUS tape
- Nb₃Sn middle coil (9.9T)
 - internal tin Nb₃Sn inner coil
 - high strength Nb₃Sn middle coil
 - bronze route Nb₃Sn outer coil
- NbTi outer most coil (6.5T)

Cooling system

- single stage GM cryocoolers
(100W@50K)x2
for radiation shield, high-T end of CL
- GM-JT refrigerator (4.2W@4.2K)
for coils, low-T end of CL

Road Map

